


# Ultimate Saliva Control

## The DryTips® Advantage Over Traditional Cotton Rolls:

- Complete coverage of parotid orifice and entire buccal mucosa with polyethylene film for cutting protection
- Facilitates work far back in the oral cavity
- Flexible; accommodates to cheek movements
- Comfortably stays in place, yet is easily removed with water spray
- Reflective coating for increased visibility and to aid in retraction



NEW! Reflective coating now available.

## Ideal for:

- Impressions
- Occlusal Adjusting
- Orthodontics
- Paint-on techniques
- Composites
- Cementation
- Fissure sealing
- Fluoride applications

## DryTips® Ordering Information:

ITEM NUMBER	SIZE	QUANTITY	COLOR
291542	Large	50	White
291543	Small	50	White
291726	Large	50	Silver
291727	Small	50	Silver

Order today through your preferred dealer.

For questions, call Microbrush® International at (866) 866-8698.

**BUY 3** Boxes of DryTips® or Reflective DryTips®

**GET 1** Box of the same for **NO CHARGE!**



Use Promo Code Q119DT on Invoice

\*To redeem these offers, send a copy of your qualifying invoice to: [freegoods@microbrush.com](mailto:freegoods@microbrush.com) or fax it to 888-771-7573. Purchases must be shown on one invoice. Offer limited to one per office and must be received no later than one month after the promotion expires. Offer good 1/1/2019-3/31/2019.

**NEW!**

# CONTACT PRO<sup>®</sup>

MATRIX SYSTEM

## Simple Technique for Outstanding Results

The reliable, consistent and easy to use system that provides solutions for Class II Restorations.

LIMITED TIME OFFER!

# BUY 1

ContactPro<sup>®</sup>  
Matrix Kit System

# GET 1

of equal or lesser value  
for **NO CHARGE!**



Featuring our exclusively engineered ContactPro<sup>®</sup> matrix ring, which provides consistent results.

Use Promo Q119KIT on Invoice

### Matrix Kit Contents

Description	Item#	ContactPro Ring	SmartView Clip	Coated Matrix Bands	Wood Wedges	Forceps	DryTips	Applicators
Expert Kit - Coated	291699	2	2	3.5mm (50 bands); 4.5mm (50 bands); 5.5mm (50 bands)	50 X-Small; 50 Small; 50 Medium	1 Band Forcep; 1 Ring Forcep	20 Small; 20 Large	15 Regular Size (2.0mm)
Primary Kit - Coated	291701	2	1	3.5mm (25 bands); 4.5mm (25 bands); 5.5mm (25 bands)	25 X-Small; 25 Small; 25 Medium	1 Band Forcep; 1 Ring Forcep	-	-
Intro Kit - Coated	291703	2	-	3.5mm (10 bands); 4.5mm (10 bands); 5.5mm (10 bands)	10 X-Small; 10 Small; 10 Medium; 10 Large	1 Ring Forcep	-	-

Order today through your preferred dealer.

For questions, call Microbrush<sup>®</sup> International at (866) 866-8698.

\*To redeem these offers, send a copy of your qualifying invoice to: [freegoods@microbrush.com](mailto:freegoods@microbrush.com) or fax it to 888-771-7573. Purchases must be shown on one invoice. Offer limited to one per office and must be received no later than one month after the promotion expires. Offer good 1/1/2019-3/31/2019.